

IN VITRO FERTILIZATION



LAMSYSTEMS

IVE
WORKSTATIONS

Specialty designed for Assisted Reproductive Technology

2020

IVF

WORKSTATIONS ARE SPECIALLY DESIGNED FOR LABORATORIES OF IN VITRO FERTILIZATION.

Highly pure air environment in the cabinet chamber provides solid protection of the technological process and minimizes the risk of contamination when working with oocytes and embryos and doing research in ART.



FEATURES


CLEAN AIRFLOW

- Consistent, turbulence free airflow with HEPA H14 filtration according to the European Standard EN 1882-1 provides air cleanliness efficiency of 99,995% for particles sized 0,3 µm
- Carbon filter as part of the standard configuration provides air purification from VOC.
- Radial low-noise EC fans decrease noise and vibration level for comfortable working conditions



COLOUR TOUCHSCREEN

- Selection and adjustment of the airflow settings (pre-operation, operation, economy, clean maintenance modes)
- Automatic switch-on time setting for the cabinet's pre-readiness to work on the specified days of the week
- Tabletop heating temperature setting (in case of two heating zones, the temperature is individually set for each one)



WORKING AREA

- L-shaped heated zone/-s
- Different location of heated zones
- Smooth single-piece tabletop made of stainless steel has a visually marked heating zone to ensure the samples are always warm



INDEPENDENT HEATED ZONES

Each workplace is featured with 7 independent heated zones with individual sensors and heating elements that ensure accurate temperature adjustment as well as fast and even heating



LIGHT AND HEATED OBJECT PLATE

- An integrated LED light with heated object plate is featured with independent temperature adjustment
- Light angle regulation
- Lighting dimmer is located on the tabletop
- LED light spectrum in the work chamber excludes UV radiation and its harmful impact

IVF

WITH LAMINAR AIRFLOW



PRODUCT PROTECTION

Highly pure air environment is created in the work chamber to protect valuable samples during ART.

- All brands of microscopes can be installed
- Additional options: gas supply and humidification system, 21.5" monitor, antivibration table



Built-in system for premixed gas supply and humidification



Cut-out for the microscope (shape, size and location as per customer's requirements)





Removable armrest



Light angle adjustment



AIR FLOW SCHEME

VARIATIONS OF HEATED SURFACE LOCATION AND ARRANGEMENT


Mirrored positions are also available.

TECHNICAL SPECIFICATION

Air cleanliness class in the work chamber for suspended particle (aerosol) concentration according to ISO 14644-1	Class 5
Class of the HEPA filters according to EN 1822 -1	H14
Primary efficiency of HEPA filters for particles sized 0.3 μm according to EN 1822 -1, %	99,995
Average downflow velocity in the work chamber in operating mode / 50% mode, m/s	0,40±0,03 / 0,25±0,03
Work zone illumination, lx, not less than	2000
Modes of Operation: Set-up - blowdown Maximum fan rotation for 1 min Operating mode..... Preset operating airflow velocity Clean maintenance mode..... Fan in economy mode, lighting and other functions are off Economy mode..... 50% fan capacity	
Warm-up by specified time	with timer and calendar

HEATED SURFACES

Preset temperature range, °C	+35 to +45
Accuracy of preset temperature maintenance, °C	±0,3

MAIN PARAMETERS AND DIMENSIONS

Article	1E-D.006-12.0	1E-D.006-15.0	1E-D.006-18.0
Dimensions of the cabinet assembled with the stand (WxDxH), mm	1200x710x1915	1500x700x1915	1800x700x1915
Dimensions of the work chamber (WxDxH), mm	1130x630x660	1425x630x640	1725x630x640
Input power of the cabinet, W, not more than (exclusively of the load on the built-in outlets)	980	980	1780
Maximum acceptable load on the built-in outlets, W, not more than	1000	1000	1000

IVEF

CLASS II TYPE A2




OPERATOR, ENVIRONMENT,
PRODUCT PROTECTION



Highly pure air environment is created in the work chamber to protect valuable samples, environment, and operator during ART.


- All brands of microscopes can be installed
- Curtain for microscope
- Three-section tabletop. Center section is fixed. Side sections are featured with gas springs for convenient disinfection
- Additional options: armrest, gas supply and humidification system, 21.5" monitor



Microscope cutout is equipped with a flap ensuring cleanliness of the work zone in clean maintenance mode



Special curtain covers the microscope cutout maintaining protection efficiency of the cabinet in operating mode



Power-driven sliding front sash



Rear wall can be featured with a built-in 21.5" monitor (additional option)



ENERGY EFFICIENCY



EC FANS

Class II microbiological safety cabinets SAVVY^{SL} are equipped with centrifugal, energy-efficient and low-noise EC FANS that significantly decrease operating costs as well as reduce the level of acoustic noise and vibration ensuring comfortable work of the personnel.

ADVANTAGES OF EC FANS:

- Monitoring and accurate adjustment of operating modes with the microprocessor control system
- Low energy consumption
- Low heat emission
- Low noise level
- No vibration
- Extended operating life

LOW POWER CONSUMPTION OF THE CABINET

0.112 kW

COMPARISON

	Input Power kW	Power Consumption per Year kWh ^[2]	CO ₂ Emissions t/year ^[3]	SAVING	CO ₂ REDUCTION
LAMSYSTEMS	0,112 ^[1]	233,0	0,123	30%	30%
Alternative*	0,160**	332,8	0,175		



* Equipment with equivalent technical characteristics produced by a known manufacturer was taken for comparison.

** Information was taken from the official advertising materials of the manufacturer.

[1] – The measurements are taken at operating mode whereat the fans and the work chamber lighting are on; the load on the built-in electric sockets is excluded.


[2] – 8 hours per day, 5 days, 52 weeks.

[3] – Each kWh of energy produced corresponds to 0.527 kg of CO₂ emission (source: <https://www.carbonindependent.org/15.html>)

LOW ENERGY CONSUMPTION PROVIDES FOR LOW HEAT EMISSION REDUCING NECESSITY IN ROOM AIR CONDITIONING AND, THEREFORE, ITS COST.

HEAT EMISSION READINGS&CALCULATIONS

Microbiological safety cabinets generate heat that may cause room temperature increase and air humidity decrease leading to operator's discomfort, loss of efficiency, fatigue, skin irritation and itching. Low energy consumption provides for low heat emission reducing necessity in room air conditioning and, therefore, its cost.





Measuring the Temperature of the Air Inflow



Measuring the Temperature of the Air Outflow




Measuring the Temperature of the Front Panel Heated Area



HEAT EMISSION OF THE CABINET:

- with fans and lighting on **45W**
- with fans on and lighting off **25 W**
- with lighting on and fans off **20 W**

COMPARISON



- Operating desktop laser printer..... 215 W
- Working person 180 W
- Resting person 100 W
- Operating monitor (19") 80 W
- IN THE OPERATING MODE*** **45 W**

* Equipment with equivalent technical characteristics produced by a known manufacturer was taken for comparison.

** Information was taken from the official advertising materials of the manufacturer.

* The measurements are taken at operating mode whereat the fans and the work chamber lighting are on; the load on the built-in electric sockets is excluded.




LOW NOISE LEVEL 47dBA

The level of noise in real operating environment depends on the dimensions of the operating site, on the cabinet's location as well as on the total background noise and may vary by 3-4 dB(A).

SIMPLE AND CONVENIENT SETTING OF AIRFLOW VELOCITY

The system ensures separate control of the inflow and downflow velocities as well as automatically maintains the air balance. There is no need to adjust the air balance manually and, therefore, the maintenance time in case of qualification, filter replacement or periodic verifications is significantly reduced.

HIGHLY ACCURATE MAINTENANCE OF PRESET AIRFLOW VELOCITY at any level of filter clogging and in case of changing ambient conditions (humidity, temperature, pressure).



FRONT SASH STOP

To ensure maximum safety, the sliding front sash automatically stops in two cases:

POWER DRIVEN FRONT SASH

Five preset positions of the front sash:



1. OPERATING POSITION in Main Operation mode.
2. UP POSITION for work chamber loading/unloading.
3. STOP BEFORE CLOSING for safety.
4. CLOSED POSITION in clean maintenance mode or for work chamber UV irradiation.
5. DOWN POSITION for disinfection of the upper part of the front sash.

BEFORE CLOSING

Automatic stop of the front sash at 40 mm height from the surface prior to complete closing of the front opening prevents any injury of the operator's hands.

AHEAD OF OBSTACLE

The moving front sash automatically stops in case of contact with an object or a hand of the operator.

AIR FLOW SCHEME

VARIATIONS OF HEATED SURFACE LOCATION AND ARRANGEMENT


TECHNICAL SPECIFICATION

Air cleanliness class in the work chamber for suspended particle (aerosol) concentration according to ISO 14644-1 Class 5

Class of the cabinet according to EN 12469, NSF/ANSI 49 II

Type of the cabinet according to NSF/ANSI 49 A2

Class of the HEPA filters according to EN 1822-1 H14

Primary efficiency of HEPA filters for particles sized 0.3 μm according to EN 1822-1, % 99,995

Average downflow velocity in the work chamber in operating mode / 50% mode, m/s 0,47±0,03 / 0,20±0,02

Work zone illumination, lx, not less than 2000

Air recirculation, % 70

Modes of Operation:

Set-up - blowdown Maximum fan rotation for 1 min

Operating mode Specified operating airflow velocity

Clean maintenance mode Fan in economy mode, lighting and other functions are off

Economy mode 50% fan capacity

Warm-up by specified time preset with timer

HEATED SURFACES

Preset temperature range, °C +35 to +45

Accuracy of preset temperature maintenance, °C ±0,3

MAIN PARAMETERS AND DIMENSIONS

Article	1E-B.006-12.0	1E-B.006-15.0
Dimensions of the cabinet assembled with the stand (WxDxH), mm	1200x810x2095	1500x810x2095
Dimensions of the work chamber (WxDxH), mm	1105x610x700	1405x610x700
Input power of the cabinet, W, not more than (exclusively of the load on the built-in outlets)	1150	1150
Maximum acceptable load on the built-in outlets, W, not more than	1000	1000
Power of the UV lamp, W	30	30

LAMSYSTEMS

LAMSYSTEMS is the largest manufacturer of high-tech equipment for personnel, product, and environment protection from biologically hazardous agents.

Company supports full cycle of production including design, manufacturing, sales, and maintenance services of specialized equipment such as microbiological safety cabinets, laminar flow benches, fume hoods, clean zones, etc. All produced equipment has received quality certificates.



www.lamsys.com



LAMSYSTEMS

LAMSYSTEMS GmbH
Magdeburger Str. 3,
14641 Wustermark bei Berlin,
Germany
Tel.: +49 (0) 30 2555 9888
info@lamsys-euro.com

BMT EXCLUSIVE DEALER

Published in 2020.

The manufacturer reserves the right to change specifications and designs in further technical development of the equipment